






INSIGHT™

MRI Overlay & Coil Holders

Enabling MR Imaging in the Same Position as Treatment



MR Imaging & CT Imaging



Through Treatment



positioning
patients for life.®

www.Qfix.com

INSIGHT™

INSIGHT™

INSIGHT™ Overlay and Coil Holders are designed to enable MR imaging in the same position as treatment. The Insight MR Overlay thin design positions the patients closer to the scanner coils for improved MR imaging. At 6.1 kgs (13.5 lbs) the lightweight design makes it easy to remove and install for variety of imaging setups. The Varian Exact® Indexing notches at 70 millimeters provides more flexibility for indexing positioning accessories.


INSIGHT™ MR OVERLAY

INSIGHT™ MR Overlay mounts to compatible Siemens® MAGNETOM™ Aera™, Skyra™, Vida™ and Sola™ table tops allowing the use of MR compatible patient positioning devices during MR scanning for radiotherapy. The INSIGHT™ MR Overlay is compatible with Portrait™ head only S-frame thermoplastic masks.


INSIGHT™ MR Coil Holders

INSIGHT™ MR Coil Holders are intended to position body and head coils allowing the patient to undergo imaging in the same position as that of the treatment position. Coil holders are available for Siemens® MAGNETOM™ Aera™, Skyra™, Vida™ and Sola™ scanners using the Flex 4, Body 6, and Body 18 coils.

Please contact Qfix® for more details on compatibility with Siemens® MR Systems.



Accessories:



Portrait™ Head Only S-Frame Masks are compatible with the INSIGHT™ MR Overlay. *See Qfix.com for more types and descriptions.*

RT-4700-01	INSIGHT MR Overlay
RT-4700-02	INSIGHT MR Head Coil Holders
RT-4700-03	INSIGHT MR Body Coil Holders
RT-4700-10	INSIGHT Complete Package (includes MR Overlay, MR Head & Body Coils)

2007559_EN_D

© 2018. All rights reserved. Qfix is a registered trademark of Qfix Systems, LLC.

This brochure contains information about products which may or may not be available in particular countries. Each country has specific laws and regulations governing the commercialization of medical devices and the communication of information regarding medical devices in printed and digital media. Products appearing in this brochure or on our website may or may not have received approval, clearance, or marketing authorization by a governmental regulatory body in any particular country, or may have received approval, clearance, or marketing authorization for different indications and restrictions in different countries, or may be for investigational use only. A product's appearance in this brochure should not be construed as a solicitation or promotion for said product, nor any indication which is not authorized by the laws and regulations of the country where the reader resides.



Qfix
Avondale, PA
+1 484.720.6053
sales@Qfix.com
www.Qfix.com