




Elevate your Workflow



SRS



Pelvis



Intracranial,
Head & Neck



Abdomen




Craniospinal



Lung



Breast





SBRT



Quickly Adapts to a Variety of Clinical Treatments

Alta is an all-in-one solution that quickly adapts to a wide range of setups, reducing setup times.

- Combines functionality of multiple devices
- Dependable and reproducible
- Fast and efficient setups
- Can be extended beyond end of couch top for improved head and neck dosimetry
- Lightweight – 5.4 kg (12 lbs)
- CT, MRI, and treatment compatible



InfinityEdge™ & InfinityLock™

InfinityLock accessories tightly fasten to the device and can be infinitely indexed along the InfinityEdge of the device.



Alta™ ArmShuttle™ Elite

Quickly and securely attaches to provide comfortable, reproducible positioning of the arms out of the treatment field.





Elevate Your Workflow



Scan to View
Alta Multipurpose
Device Video

Alta Standard Package

RT-4800-10




Knee & Foot Supports

- Ergonomically designed for patient comfort while maintaining reproducible positioning
- Optional clip-on leg and foot straps for added immobilization
- Highly durable and easily cleanable



Thoracic & Pelvis

Alta allows positioning and immobilization of patients for thoracic and pelvic applications through the use of thorax and pelvis thermoplastics and other pelvic positioning products.



RT-4800-GL
Alta Groin Lock



RT-4800-H45
Alta Opti-Handle for Pelvis, 45cm, set of 2



RT-4800-H30
Alta Opti-Handle for Pelvis, 30cm, set of 2

Stereotactic Radiosurgery (SRS)


Alta can easily be used as a stereotactic radiosurgery treatment solution through the quick attachment of the SRS Adapter.



RT-4600-03-MRI
MR SRS Adapter



SRS Headrests Included



Stereotactic Body Radiation Therapy (SBRT)

The low electron generating properties of the Alta Multipurpose Device make it the ideal solution for a range of radiotherapy indications from palliative treatment to stereotactic body radiation therapy (SBRT).



RT-4800-B1
InfinityLock SBRT Compression Bridge Small (min: 30 cm & max: 38 cm)

RT-4800-B2
InfinityLock SBRT Compression Bridge Large (min: 39 cm & max: 48 cm)

RT-4558KV4-01
Low Profile Compression Paddle

RT-4558KV4
Compression Paddle



RT-4558CB03
ZiFix Traverse Abdominal/Thoracic Motion Control System



Qfix
Avondale, PA
+1 484.720.6053
sales@Qfix.com
www.Qfix.com

2008576_EN_D

This brochure contains information about products which may or may not be available in particular countries. Each country has specific laws and regulations governing the commercialization of medical devices and the communication of information regarding medical devices in printed and digital media. Products appearing in this brochure or on our website may or may not have received approval, clearance, or marketing authorization by a governmental regulatory body in any particular country, or may have received approval, clearance, or marketing authorization for different indications and restrictions in different countries, or may be for investigational use only. A product's appearance in this brochure should not be construed as a solicitation or promotion for said product, nor any indication which is not authorized by the laws and regulations of the country where the reader resides. © 2022 All rights reserved. Qfix is a registered trademark of Qfix. Alta, ArmShuttle, FlexLock, InfinityEdge, InfinityLock, and ZiFix Traverse are trademarks of Qfix.